Prior to use, carefully read the instructions for use.

**EN**

## GC Fuji PLUS

RADIOPAQUE REINFORCED GLASS IONOMER LUTING CEMENT

For use only by a dental professional in the recommended indications.

#### RECOMMENDED INDICATIONS

- Cementation of metal-based inlays, onlays, crowns and bridges
- Cementation of resin inlays, onlays, crowns and bridges
- Cementation of all ceramic inlays
- Cementation of high strength (zirconia based) all ceramic crowns and bridges
- Cementation of metal, ceramic and fiber posts

#### CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

#### DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	2,0/1,0
Mixing time (sec.)	20"
Working time (min., sec.) (at 23 <span> </span> °C) (from start of mixing)	2' 30"
Final finishing commencing time (after seating the restoration)	4' 30"

#### 1. TOOTH PREPARATION

- Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
  - Clean the prepared tooth with pumice and water.

Note:

Use GC Fuji PLUS CONDITIONER on the prepared tooth surface for 20 seconds for optimal adhesion.
  - Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE. Prepared surfaces should appear moist (glistening).
- RESTORATION PREPARATION

Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.

- POWDER AND LIQUID DISPENSING
  - The standard powder to liquid ratio is 2.0g / 1.0g. (1 level large scoop of powder to 3 drops of liquid, or 1 level small scoop of powder to 1 drop of liquid).
  - For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.
  - Hold the liquid bottle vertically and squeeze gently. d) Close bottles immediately after use.

- MIXING

Dispense powder and liquid onto the pad. Using the plastic spatula, add all the powder to the liquid. Mix rapidly for 20 seconds.

Note :
When mixing larger amounts, divide the powder into 2 equal parts. Mix the first portion with all of the liquid for 15 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-25 seconds (total 30-40 seconds).

- CEMENTATION
  - Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 30 seconds from start of mixing at 23°C (73,4°F). Higher temperature will shorten working time.

Note:

For more complex procedures requiring longer working time, it is recommended to mix powder and liquid on a cold mixing slab and or use GC Fuji PLUS EWT.
  - Maintain moderate pressure.
  - Start removing excess cement when the excess cement feels rubbery.
  - Finishing can be started 4 minutes 30 seconds after seating the restoration.

**STORAGE**  
Store in a cool and dark place (4-25°C) (39.2-77.0°F). (Shelf life: powder 3 years, liquid 2 years from date of manufacture)


- SHADES**  
Yellow  
A3 (Vita® shade)  
Translucent  
Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.
- PACKAGES**  
1. 1-1 package : 15g powder, 8g (7mL) liquid, GC Fuji PLUS CONDITIONER 7g (6.5mL), a double-ended powder scoop, a plastic spatula, mixing pad (#23)  
2. Powder 15g with scoop  
3. Liquid 8g (7mL)  
4. GC Fuji PLUS CONDITIONER 7g (6.5mL)

**CAUTION**

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with components of other glass ionomer cements.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
- Do not use GC Fuji PLUS in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS from setting or bonding properly.
- This product is not indicated for filling or core build-up.

Last revised : 04/2008

--

--

MANUFACTURED BY  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED BY  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

 EU: GC EUROPE N.V.  
Research-Park, Interleuvenlaan 13, B-3001 Leuven, Belgium  
TEL: +32-16-39.80.50

GC AMERICA INC.  
3737 West 127th Street, Alsip, IL 60803 U.S.A.  
TEL: +1-708-597-0900

GC ASIA DENTAL PTE. LTD.  
19 Loyang Way, #06-27 Singapore 508724  
TEL: +65 6545 7588

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