Prior to use, carefully read the instructions for use.

EN

## G-CEM CAPSULE

SELF ADHESIVE LUTING CEMENT IN CAPSULES

For use only by a dental professional in the recommended indications.

### RECOMMENDED INDICATIONS

- Cementation of all types of all ceramic, resin and metal-based inlays, onlays, crowns and bridges (max. 2 abutments)
- Cementation of metal, ceramic, fiber posts, and cast post and cores.

### CONTRAINDICATIONS

- Pulp capping.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

### DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	0,35/ 0,15
Standard Mixing Time (sec.) (4000 rpm at 23°C)	10"
Working Time (min., sec.) (at 23°C)	2'30"

Test conditions :Relative humidity (50+/-10%).

### A. CEMENTATION TECHNIQUE for ceramics, indirect composites and metal based restorations

#### 1. TOOTH PREPARATION

- Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
- Clean the prepared tooth with pumice and water.
- Remove the tooth with water. Dry by gently blowing with oil free air. Prepared surfaces should appear dry.

NOTE:

- Do not use substances like H<sub>2</sub>O<sub>2</sub> (hydrogen peroxide) and/or disinfectants as these may lower the bond strength.
- Special care should also be taken to remove so called imaging powders very thoroughly with pumice and water.

#### 2. RESTORATION PREPARATION

Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.

#### 3. CAPSULE ACTIVATION AND MIXING

- After activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1).
- To activate the capsule, push the plunger until it is flush with the main body (Fig. 2).
- Immediately place the capsule into a GC CAPSULE APPLIER and click the cover once (Fig. 3). The capsule is now activated.

- The capsule should be activated just before mixing and used immediately.
- Immediately remove the capsule and set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4.000RPM) (Fig. 4). More mixing on mixing times can be found in following chart.

#### Mixing time chart

Mixer Speed	Mixing Time
Standard mixer (approx 4.000rpm)	10 seconds
High-speed mixer (over 4.500rpm)	

Rotational mixer  
\*No centrifuge option

Note: Refer to the respective manufacturer's instructions for use.

#### 4. DISPENSING

- Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.
- Make two clicks to prime the capsule (Fig. 5). To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.

- Extrude directly into the restoration (Fig. 6).
- Note:  
Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation procedure.

- To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.

#### 5. CEMENTATION

- Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 30 seconds from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time. Refrigeration of G-CEM CAPSULE is recommended to extend the working time.

- Maintain moderate pressure.
- EXCESS CEMENT REMOVAL

- To remove excess cement, light cure using a visible light curing unit (Halogen, GC G-Light, LED) for 2-4 seconds (Fig. 7) or allow the cement to self-cure until the excess cement feels rubbery.

- Light cure each surface for 10 seconds (single surface from occlusal:10 seconds, any other surfaces:10 seconds each) or allow the cement to self cure for 4 minutes after seating of the restoration.

#### NOTE:

- G-CEM is a dual cure cement. To insure proper curing involving cementation of the following restorations, G-CEM must be used in self cure mode for 4 minutes.  
-Metal, metal coping restorations  
-Metal-free dark and opaque restorations  
-Translucent metal-free restorations greater than 2mm in thickness

- Finish and polish margins with appropriate instruments.

#### B. CEMENTATION TECHNIQUE for metal, fiberglass and ceramic posts

##### 1. TOOTH PREPARATION

Prepare the post space according to manufacturer's instructions. The post space should be prepared and dried thoroughly using paper points.

- Make sure that the post is pretreated and handled according to the manufacturer's instructions.

##### 3. CAPSULE ACTIVATION AND MIXING

- After activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1).
- To activate the capsule, push the plunger until it is flush with the main body (Fig. 2).
- Immediately place the capsule into a GC CAPSULE APPLIER and click the cover once (Fig. 3). The capsule is now activated.

- The capsule should be activated just before mixing and used immediately.
- Immediately remove the capsule and set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4.000RPM) (Fig. 4). More mixing on mixing times can be found in following chart.

#### Mixing time chart

Mixer Speed	Mixing Time
Standard mixer (approx 4.000rpm)	10 seconds
High-speed mixer (over 4.500rpm)	

Rotational mixer  
\*No centrifuge option

Note: Refer to the respective manufacturer's instructions for use.

#### 4. DISPENSING

- Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.
- Make two clicks to prime the capsule (Fig. 5). To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.

- Extrude directly into the restoration (Fig. 6).
- Note:  
Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation procedure.

- To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.

#### 5. CEMENTATION

- Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 30 seconds from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time. Refrigeration of G-CEM CAPSULE is recommended to extend the working time.

- Maintain moderate pressure.
- EXCESS CEMENT REMOVAL

- To remove excess cement, light cure using a visible light curing unit (Halogen, GC G-Light, LED) for 2-4 seconds (Fig. 7) or allow the cement to self-cure until the excess cement feels rubbery.

- Light cure each surface for 10 seconds (single surface from occlusal:10 seconds, any other surfaces:10 seconds each) or allow the cement to self cure for 4 minutes after seating of the restoration.

#### NOTE:

- G-CEM is a dual cure cement. To insure proper curing involving cementation of the following restorations, G-CEM must be used in self cure mode for 4 minutes.  
-Metal, metal coping restorations  
-Metal-free dark and opaque restorations  
-Translucent metal-free restorations greater than 2mm in thickness

- Finish and polish margins with appropriate instruments.

#### B. CEMENTATION TECHNIQUE for metal, fiberglass and ceramic posts

##### 1. TOOTH PREPARATION

Prepare the post space according to manufacturer's instructions. The post space should be prepared and dried thoroughly using paper points.

- Make sure that the post is pretreated and handled according to the manufacturer's instructions.

##### 3. CAPSULE ACTIVATION AND MIXING

- After activation, shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 1).
- To activate the capsule, push the plunger until it is flush with the main body (Fig. 2).
- Immediately place the capsule into a GC CAPSULE APPLIER and click the cover once (Fig. 3). The capsule is now activated.

- The capsule should be activated just before mixing and used immediately.
- Immediately remove the capsule and set it into a mixer (or an amalgamator) and mix for 10 seconds (+/- 4.000RPM) (Fig. 4). More mixing on mixing times can be found in following chart.

#### Mixing time chart

Mixer Speed	Mixing Time
Standard mixer (approx 4.000rpm)	10 seconds
High-speed mixer (over 4.500rpm)	

Rotational mixer  
\*No centrifuge option

Note: Refer to the respective manufacturer's instructions for use.

#### 4. DISPENSING

- Immediately remove the mixed capsule from the mixer and load it into the GC CAPSULE APPLIER.
- Make two clicks to prime the capsule (Fig. 5). To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.

- Extrude directly into the restoration (Fig. 6).
- Note:  
Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation procedure.

- To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.

#### 5. CEMENTATION

- Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 30 seconds from start of mixing at 23°C (73.4°F). Higher temperatures will shorten working time. Refrigeration of G-CEM CAPSULE is recommended to extend the working time.

- Maintain moderate pressure.
- EXCESS CEMENT REMOVAL

- To remove excess cement, light cure using a visible light curing unit (Halogen, GC G-Light, LED) for 2-4 seconds (Fig. 7) or allow the cement to self-cure until the excess cement feels rubbery.

- Light cure each surface for 10 seconds (single surface from occlusal:10 seconds, any other surfaces:10 seconds each) or allow the cement to self cure for 4 minutes after seating of the post (Fig. 11).

#### STORAGE

Store the capsules in the original aluminium foil in a cool and dark place (4-25°C (39-77.0°F)). (Shelf life: 2 years from date of manufacture)

#### SHADES

A2 (Vita® shade), A03 (opaque), B01 (opaque) and Translucent  
\*Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

#### PACKAGING

- Single shade package : 50 capsules.
- Assorted package : 50 capsules (20 each of A2 and Translucent, 5 each of A03 and B01).

- Average contents per capsule : 0.35g powder and 0.15g (0.14mL) liquid.
- Minimum net volume of mixed cement per capsule : 0.17mL.

- GC CAPSULE APPLIER (1 piece) is available separately.
- GC CAPSULE Elongation Tip: 50 pieces

#### CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, use a rubber dam and/or cocoa butter can be used to isolate the operation field from the skin or oral tissue.

- In case of contact with eyes, flush immediately with water and seek medical attention.
- If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mark disappears, usually in 1-2 weeks.

- Do not use G-CEM in combination with eugenol containing materials as eugenol may hinder G-CEM from setting or bonding properly.
- Use a protective light shield or similar protective eye wear during light curing. This product is not indicated for filling or core build-up.
- When wiping elongation tip with a cotton soaked with alcohol, prevent the alcohol from entering the inside of the tip. Elongation tip is for single use only and cannot be autoclaved or chemically.

Last revised : 02/2011

0086

IGI

MANUFACTURED BY  
GC CORPORATION  
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

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